

# Lovibond® Water Testing

## Tintometer® Group



### Step 1

Please store the update file "bod.hex" on a SD-card or USB-stick. Make sure that the file appears at first in your root directory to enable your instrument recognizing the file.

### Step 2

Turn on your instrument in mains operation and connect it to the data carrier.

### Step 3

Choose "update" in menu options and start the update with the right arrow key. Confirm the security question with OK and wait until the boot loader finished.

### Step 4

The device restarts automatically and shows the main menu.

### Step 5

After all please check the revision no. in menu option "update".

## Instruction Software-Update

# Lovibond® Water Testing

## Tintometer® Group



## Changelog

### New regarding MD600/MultiDirect:

- New Methods Bromine PP,  $\text{ClO}_2$  PP,  $\text{H}_2\text{O}_2$  LR and  $\text{H}_2\text{O}_2$  HR
- New way to handle M265 Nitrate tube test

---

### New regarding PM620/PoolDirect:

- New Method  $\text{H}_2\text{O}_2$  HR