

# Lovibond® Water Testing

## Tintometer® Group



### Step 1

Please open the instruction „HexLoad\_de\_gb.pdf“ in the „HexLoad“ dir.  
Print the instruction and have it at hand when updating.

### Step 2

To update your Photometer or Turbidimeter  
open the program HEXLoad.

### Step 3

Additional every Product has its own firmware. You will find it in the folder  
„Firmware“ fitting to your product.

### Step 4

Update your device as described in the instruction „HexLoad\_de\_gb.pdf“.

### Step 5

Download the current version of the instruction for your product from our website.

## Instruction Software-Update